



FT-707

DIGITAL THERMOMETER

- Digital Thermometer
- Clinical accuracy
- Water proof
- Auto shut-off
- Fever Alert
- Break resistance
- Beeps when complete
- Memory recall
- For use at multiple sites
- BPA Free, Latex Free
- No mercury or glass
- Memory of last measured test
- With protective cover
- Medical product

- Size: approx. 1.7 x 12.4 x 1.0 cm
- Weight: approx. 11g, including battery
- Range: 90.0°F - 109.9°F
- Accuracy: $\pm 0.2^{\circ}\text{F}$
- Battery: 1 x 1.55V LR41

- EAN no.: 747 3800 70700 5
- Pack: 1 (10xFT-707)/ Shipping carton:
30x1

