

CERTEZA[®]
...your good health



NB-607

NEBULIZER COMPRESSOR SYSTEM

- Compressed air nebulization
- For treatment of upper and lower airways, colds, asthma and respiratory diseases
- Premium carrying bag
- Optimal respirable dose delivery
- Quick effective treatment
- Support for wide voltage input
- Durable quality product
- Non-slip body material
- Built-in medicine cup holder
- Angled-mouth piece
- Long tube
- Medical product
- Size: approx. 14.5 x 19.0 x 9.0 cm
- Weight: approx. 1.4kg (machine only)
- Power output: 57 Watt
- Nebulizer chamber capacity: 10ml max
- Nebulization rate: 0.27 ml/min
- Particle size: 3 μ m (MMD)
- EAN no.: 747 6322 60700 6
- Pack: 6/shipping carton

