



**CERTEZA®**  
...your good health

**BM-400**

## ARM BLOOD PRESSURE MONITOR

Blood pressure and pulse measurement on arm  
Measurement during inflation  
New attractive design  
For arm circumference of 22-40 cm  
USB-C external power source port  
Movement error indicator  
Large LCD display  
Hypertension risk indicator  
2 x 120 memory recall function  
Irregular Heartbeat Indicator  
Average of last 3 readings  
One touch automatic operation  
Date and Time Function  
Automatic switch off  
Signal in case of application error  
Low battery indicator  
Storage pouch  
Medical product

For arm circumference: 22 – 40 cm  
Size: approx. 12.6 x 8.5 x 5.7 cm  
Weight: approx. 148g, excluding batteries  
Measuring method: Oscillometric method  
Measuring range:  
– Blood pressure: 60 mmHg – 230mmHg  
– Pulse: 40 – 170 beats/min  
Batteries: 3 x AA (included)  
USB-C adapter : 5V 1A (optional)  
Pouch for storage

EAN no.: 747 8276 40000 9  
Pack: 20 pcs / Shipping carton

