

CERTEZA[®]
...your good health



98.6

99.5

102

FT-712

INFRARED FOREHEAD THERMOMETER

No skin contact
Reliable results in one seconds
For forehead and objects
Safe and hygienic for babies
3 Colour LCD fever indicator
30 memory spaces
Selectable °C or °F
Date and time display
Large LCD with backlight
Easy to clean
2xAAA batteries included
Medical product
Professional accuracy

Size: approx. 16.8 x 3.79 x 4.73 cm

Weight: approx. 97g, including battery

Range:

- Forehead 34°C to 43°C (93.2°F to 109.4°F)

- Object 0°C to 100°C (32°F to 212°F)

Accuracy:

- Forehead mode $\pm 0.4^{\circ}\text{F}$ or $\pm 0.2^{\circ}\text{C}$ (conditional)

- Object mode $\pm 4\%$ or $\pm 2^{\circ}\text{C}$ ($\pm 4^{\circ}\text{F}$)

Battery: 2 x AAA (included)

EAN no.: 747 6228 71200 2

Pack: 90 pcs / shipping carton

