



**CERTEZA**<sup>®</sup>  
...your good health



**PS-812**

## DIGITAL PLASTIC SCALE

Full ABS plastic body  
TARE function for mother and child weighing  
High loading capacity: 200KG or 440 lb  
Large Platform: 32x33cm  
Four precision sensors  
Inverse backlight display  
Large LCD display: 75 x 32mm  
Graduation:  
    50g (<100kg)  
    100g (>100kg)  
Automatic switch-on and switch-off  
Slim design  
Overload indicator  
Low battery indicator  
Included 2xAAA Batteries  
Measuring units: Pound (lb) or Kilograms (kg)

Size: approx. 32.0 x 33.0 cm  
Weight: approx. 2.4kg (with batteries)  
Range: Upto 200kg  
Accuracy:  $\pm 0.2$  lbs or  $\pm 0.1$  kg  
Digit size: 32mm  
Battery: 2 x AAA

EAN no.:  
    Black: 747 3770 81200 6  
    White: 747 3770 81290 7  
Pack: 4/Shipping carton:

