

CERTEZA®
...your good health



PO-908

FINGER PULSE OXIMETER

Measurement of arterial oxygen saturation (SpO2) and heart rate (Pulse)

Bright TFT display

Pulse wave display

Data review table

Rotation of display in all directions

High quality; High accuracy sensor

Interference resistance against ambient light

Accurate performance at low perfusion

Low power consumption

For use at home and on the go

For persons with: Respiratory disorders, Heart failure.

Also for sportsman going at high altitude

Audio and video alarm function

Alarm and beep adjusted in settings menu

Automatic switch off

Low battery indication

Medical product

Incl. 2 x 1.5V AAA batteries, carrying pouch and lanyard

Size: approx. 6.2(L) x 3.5(W) x 3.4(H) cm

Weight: approx. 50g, with batteries

Measuring range:

- SpO2: 0% - 100%

- Pulse: 25 - 250 beats/min

Batteries: 2 x AAA

EAN no.: 747 7869 90800 7

Pack: 100/Shipping carton

