

CERTEZA[®]
...your good health



NB-605

NEBULIZER COMPRESSOR SYSTEM

Compressed air nebulization

For treatment of upper and lower airways, colds, asthma and respiratory diseases

Optimal respirable dose delivery

Quick effective treatment

Support for wide voltage input

Durable quality product

Non-slip body material

Built-in medicine cup holder

Angled-mouth piece

Power cord with fuse

Long tube

Medical product

Size: approx. 10.0 x 18.5 x 18.5 cm

Weight: approx. 1.4kg (machine only)

Power output: 57 Watt

Nebulizer chamber capacity: 8ml max

Nebulization rate: 0.27 ml/min

Particle size: 4µm (MMD) | MMAD 1.3µm

EAN no.: 747 6322 60500 2

Pack: 6/shipping carton

