

CERTEZA[®]
...your good health



GL-100

BLOOD GLUCOSE MONITOR

Code Free

Enzyme: Glucose dehydrogenase

Testing time 5 seconds

Only a tiny blood sample of 0.6ul needed

Shows blood hematocrit value

Hematocrit range : 0 - 70%

Alarm alert of testing schedule

360 test memory

7, 14, 28 and 90 days averaging

Measurement range: 9 - 600 mg/dL

Strip ejector function in meter

Long strip expiry of 24-months

Conforms to new ISO standard 15197:2013

Economical test strips

Large strip size

Custom designed lancing device

Blood sample: capillary whole blood

For in Vitro diagnostic use

For self-testing and professional use

Medical product

Size: 9.0 x 6.1 x 1.8 cm

Weight: 51.6 g, without batteries

Batteries: 1 3.0V CR2032 (included)

Battery life: 1000 tests

EAN no.: 747 2476 10000 0

Pack: 25 pcs / Shipping carton

