

CERTEZA[®]
...your good health



FT-715

FOREHEAD INFRARED THERMOMETER (GUN-TYPE)

No skin contact - hygienic and safe
Reliable results within a second
IP22 rated - Protection against dripping water and dust
Works in low temperatures (5°C and above). Suitable for school, restaurants, hospitals etc and other outdoor places
Distance sensor (blue light) for accurate measurement
Included 9 Volts large battery for long lasting usage
99 sets of memory space
3 colour LCD display with fever alert beep
Selectable °C or °F scale
Body, Object and Room temperature measurement
Sensor cover
Pure copper probe for higher accuracy
Medical device

Size: approx. 13.2(L) x 4.2(W) x 17.3(H) cm
Gross Weight: approx. 192g (With battery)

Range:

- Forehead 32°C to 43°C (89.6°F to 109.4°F)
- Object 0°C to 100°C (32°F to 212°F)
- Room 5°C to 40°C (41°F to 104°F)

Accuracy: $\pm 0.3^{\circ}\text{C}$ (conditional)

Battery: 9V DC

EAN no.: 747 6228 71500 3

Pack: 30/Shipping carton

