

CERTEZA®
...your good health



AM-205

AIR MATTRESS WITH PUMP

- Bedsore prevention and treatment
- Alternating pressure mattress
- Silent operation pump
- Pressure adjustment knob
- Hook for hanging pump
- End flaps on mattress
- Size after inflation: 90 x 200 x 6.5 cm
- Weight: 2kg
- Max capacity: 110 kg
- Mattress colour: Beige
- Air output: >4.5 lpm
- Pressure range: 45-110 mmHg
- Cycle time: 12 min
- Material: PVC
- Power output: 7 watts
- EAN no.: 747 2697 20700 1
- Pack: 4/shipping carton

