

CERTEZA[®]
...your good health



NB-603

NEBULIZER COMPRESSOR SYSTEM

Compressed air nebulization
Adjustable flow rate using kit
Low sound level of motor
Oil-free piston type
For treatment of upper and lower airways, colds,
asthma and respiratory diseases
Optimal respirable dose delivery
Quick effective treatment
Support for wide voltage input
Durable quality product
Non-slip body material
Built-in medicine cup holder
Angled-mouth piece
Power cord with fuse
Long tube
Medical product

Size: approx. 20.5 x 10.0 x 17.0 cm
Nebulizer chamber capacity: 8ml max
Nebulization rate: 0.20 ml/min (Adjustable)
Particle size: 4 μ m (MMD) | 1.80 μ m (MMAD)
Fine Particle Fraction (%) (<4.7 μ m) : >60 %
Noise level: < 65 dBA
Electrical requirement: 220V / 50Hz
Power consumption: Below 70 watts
Max. Compressor Pressure: >30 psi
Liter Flow Range: 8 LPM

EAN no.: 747 6322 60300 8
Pack: 8 pcs / Shipping carton

