

CERTEZA[®]
...your good health



BM-307

WRIST BLOOD PRESSURE MONITOR

- Blood pressure and pulse measurement on wrist
- For wrist circumference of 13.5 - 21.5 cm
- Large LCD display
- WHO classification of hypertension indicator
- 120 memory recall function
- Irregular Heartbeat Indicator
- Average of last 3 readings
- One touch automatic operation
- 2x60 memory spaces
- Medical product
- Date and Time Function
- Automatic switch off
- Signal in case of application error
- Low battery indicator
- For wrist circumference: 13.5 - 21.5 cm
- Size: approx. 6.4 x 7.7 x 3.2 cm
- Weight: approx. 115g, excluding batteries
- Measuring method: Oscillometric method
- Measuring range:
 - - Blood pressure: 0mmHg - 300mmHg
 - - Pulse: 30 - 180 beats/min
- Batteries: 2 x AAA
- EAN no.: 747 9276 30700 8
- Pack: 24/shipping carton

