

**NEW**

Type C  
Cuff Wrap alert  
Movement indicator

**CERTEZA**<sup>®</sup>  
...your good health



**BM-405**

## ARM BLOOD PRESSURE MONITOR

Blood pressure and pulse measurement on arm

Universal cuff of 22 - 42 cm

Soft cuff material

3rd generation inflation type

USB Type C Power port

USB Type C cable included

Movement error indicator

Cuff wrap guide alert

Hypertension Indicator

Average value calculation

Memory with date and time

2 users with 99 memories each

Large LCD display

Irregular heartbeat indicator

Long-life alkaline batteries

Medical product

Automatic switch off

Signal in case of application error

Low battery indicator

Carrying bag

For arm circumference: 22 - 42 cm

Size: approx. 11.0 x 11.0 x 4.0 cm

Weight: approx. 284g, excluding batteries

Measuring method: Oscillometric method

Measuring range:

- Blood pressure: 0mmHg - 300mmHg

- Pulse: 40 - 199 beats/min

Batteries: 4 x AAA

Type C cable included

EAN no.: 747 8276 40500 4

Pack: 10/shipping carton

