

CERTEZA[®]
...your good health



MG-220

BODY MASSAGER

Replicates Shiatsu Massage
Includes 3 therapeutic attachments : Three-Point,
Soothing nodes and Deep tissue
Variable intensity function
Activates muscles
Pre and Post workout muscle recovery
Enhances blood circulation
Relieves pain and discomfort
Promotes relaxation

Size: 6.5 x 4.25 x 6.65 inches
Weight: 0.84 kg approx
Power supply: 220V
Rated power: 25W
Rotary Speed: 3000 RPM

EAN no.: 747 6747 22000 2
Pack: 20 pcs / Shipping carton



CERTEZA[®]
...your good health



MG-240

BODY MASSAGER

With Heat Therapy function
Specially designed for gentle massage of head, neck and face. And other body areas.
Includes 4 therapeutic attachments : Heat, Scalp, Facial and Spot
Two Heat settings
Two Intensity levels
Relieves pain and discomfort
Pre and Post workout muscle recovery
Reduces inflammation
Joint pain relief
Enhances blood circulation
Promotes relaxation

Size: 9.3 x 4.5 x 2.6 inches
Weight: 0.62 kg approx
Power supply: 220V
Rated power: 17W
Massager Vibration Speed:
Low 3000 RPM
High 7200 RPM

EAN no.: 747 6747 24000 0
Pack: 20 pcs / Shipping carton





MG-250

BODY MASSAGER

Tapping type massage
Includes 3 therapeutic attachments : Multipoint, Spot
and Deep tissue
Variable Intensity levels
Penetrating and soothing tapping massage
Relieves pain and discomfort
Pre and Post workout muscle recovery
Enhances blood circulation
Promotes relaxation

Size: 15.5 x 6.7 x 6.0 inches

Weight: 0.95 kg approx

Power supply: 220V

Rated power: 25W

Massager Vibration Speed:

Low 1500 RPM

High 2800 RPM

EAN no.: 747 6747 25000 9

Pack: 8 pcs / Shipping carton

