

CERTEZA®
...your good health



PLUS

BLOOD GLUCOSE MONITOR

Code Free
Enzyme: Glucose dehydrogenase
Interpretation of results
Backlight display
Testing time 5 seconds
Only a tiny blood sample of 0.6ul needed
Shows blood hematocrit value
Hematocrit range : 0 - 70%
Alarm alert of testing schedule
360 test memory
7,14, and 28 days averaging
Measurement range: 9 - 600 mg/dL
Strip ejector function in meter
Long strip expiry of 24-months
Conforms to new ISO standard 15197:2013
Economical test strips
Large strip size
Custom designed lancing device

Blood sample: capillary whole blood
For in Vitro diagnostic use
For self-testing and professional use
Medical product

Size: 7.8 x 6.4 x 1.4 cm
Weight: approx 40 g, without batteries
Batteries: 2 x 3.0V CR2032 (included)
Battery life: 1000 tests

EAN no.: 747 2476 75200 5
Pack: 25 pcs / Shipping carton

