



CERTEZA®
...your good health

SI-515

STEAM INHALER

Nebulization of liquids such as saline solution and essential oils (Eucalyptus Oil etc)
Recommended for home-management of mild cases of COVID-19 symptoms
Germ-free steam inhalation
Soft comfortable mask for targeted delivery of pure mist
Adjustable steam control
Easy 5-15 minute treatment
Fast relief
Auto shut-off
Helps alleviate congestion from cold and allergies
Add essential oil to provided pad for soothing comfort
Compatible with after-market Menthol pads for soothing relief
Relieve facial pressure and symptoms from cold, cough, allergies or flu

Size: approx. 11.0 x 11.0 x 30.0 cm
Weight: approx. 660g
Power output: 130 watts

EAN no.: 747 7846 51500 6
Pack: 24/Shipping carton

STEAM NEBULIZATION

Steam as a natural aerosol, when used with essential oils becomes therapeutic steam and provides soothing relief from Coughs, Colds, Allergies and other Respiratory Symptoms

