



BR-540

SINGLE ELECTRIC BREAST PUMP

Four Modes (Stimulation, Expression, Bionic sucking and variable frequency Expression)

Nine Levels of Vacuum for stimulation and expression

All parts that come in contact with milk are made without BPA

Run directly with Type C port. (cable included)

Easy to clean and assemble

Designed for comfortable sitting position

Compact and easy to hold

Automatic Power Off

Lightweight and portable for on-the-go use

Includes all you need to express, store and feed

Size: approx. 100 x 100 x 46 mm

Gross weight: approx. 550 g

Max Vacuum (approx.): -290 mmHg

Breast shield size: 24mm

Input : 5V, 2A

EAN no.: 747 3277 55000 5

Pack: 8 pcs /Shipping carton

