

**CERTEZA**<sup>®</sup>  
...your good health



**NB-650**

## PORTABLE MESH NEBULIZER

Vibrating MESH membrane technology  
For treatment of upper and lower airways, colds,  
asthma and respiratory diseases  
Optimal respirable dose delivery  
Suitable for children and travelling  
Super quiet function  
Type C port (cable included)  
Removable AA batteries (Not included)  
Easy to clean  
Complete kit included  
Quick effective treatment  
Durable quality product  
Non-slip body material  
Medical product

Size: approx. 6.5 x 4.6 x 11.0 cm  
Gross Weight: approx. 0.31 kg  
Nebulizer chamber capacity: 14ml max  
Nebulization rate: 0.15 - 0.9 ml/min  
Particle size: <5 $\mu$ m (MMD)  
Fine Particle Fraction (%) (<4.7 $\mu$ m) : 92 %  
Power Source: DC 1.5 x 2 AA (Not included)  
DC 5V, 1A with Type C AC Adapter  
Power consumption: Below 4 watts

EAN no.: 747 6374 65000 5  
Pack: 30 pcs / Shipping carton

