



TS-100

BLOOD GLUCOSE STRIPS

Code Free

Enzyme: Glucose dehydrogenase

Testing time 5 seconds

Only a tiny blood sample of 0.6ul needed

Shows blood hematocrit value

Hematocrit range : 0 - 70%

360 test memory

7,14, 28 and 90 days averaging

Measurement range: 9 - 600 mg/dL

Long strip expiry of 24-months

Conforms to new ISO standard 15197:2013

Economical test strips

Large strip size

Blood sample: capillary whole blood

For in Vitro diagnostic use

For self-testing and professional use

Medical product

EAN no.: 747 2487 10000 0

Pack: 100 boxes / Shipping carton

Packing size :

2 x 25 strips

1 x 50 strips

