

CERTEZA[®]
...your good health



BM-450

ARM BLOOD PRESSURE MONITOR

Blood pressure and pulse measurement on arm
Attractive design with LED display
Practical cuff storage space
Voice output of results (default-off)
Measurement during inflation
For arm circumference of 22-40 cm
USB-C external power source port
Movement error indicator
Large size LED display
Hypertension risk indicator
2 x 120 memory recall function
Irregular Heartbeat Indicator
Average of last 3 readings
One touch automatic operation
Date and Time Function
Automatic switch off
Signal in case of application error
Low battery indicator
Medical product

For arm circumference: 22 – 40 cm
Size: approx. 11.8 x 11.0 x 5.2 cm
Weight: approx. 256g, excluding batteries
Measuring method: Oscillometric method
Measuring range:
– Blood pressure: 60 mmHg – 230mmHg
– Pulse: 40 – 170 beats/min
Batteries: 3 x AA (included)
USB-C adapter : 5V 1A (optional)

EAN no.: 747 8276 45000 4
Pack: 20 pcs / Shipping carton

