

CERTEZA[®]
...your good health



FT-709

DIGITAL FLEXIBLE THERMOMETER

Digital Flexible Thermometer
New attractive design
Fast reading time of 30 seconds
Extra large display
°F or °C settings
Clinical accuracy
Water proof
Auto shut-off
Fever Alert
Break resistance
Beeps when complete
Memory recall
For use at multiple sites
BPA Free, Latex Free
No mercury or glass
Memory of last measured test
With protective cover
Medical product

Size: approx. 14.3 x 2.5 x 1.4 cm
Weight: approx. 16g, including battery
Range: 90.0°F - 111.9°F
Accuracy: $\pm 0.2^{\circ}\text{F}$ or $\pm 0.1^{\circ}\text{C}$
Battery: 1 x 1.5V LR41 or SR41, UCC 392

EAN no.: 747 3830 70900 0
Pack: 1 (10xFT-708)/ Shipping carton: 30x1

