


CERTEZA[®]
...your good health



SM-500

PORTABLE PHLEGM SUCTION MACHINE

Designed to remove thick liquid such as phlegm and blood

High flow, high vacuum and low noise operation
Step-less pressure adjustment to set any pressure level

Oil-free piston pump

Concealed storage bottle design

Overflow protection designed to prevent mucus and blood leaking into machine

Used for suction in professional settings or at home

Portable design suitable for emergency use

Power: AC 220~230V 50Hz

Max negative pressure: >0.075 MPa

Pumping flow rate: > 18 L/min

Reservoir volume: 1000 mL 1 bottle

Working system: 30 min ON. 30min OFF

Unit size: 29.3 x 19.2 x 24.5 cm

Weight (approx): 3.7kg

EAN no.: 747 7776 50000 4

Pack: 4/Shipping carton

