



CERTEZA®
...your good health

GL-110

BLOOD GLUCOSE MONITOR

Code Free

Enzyme: Glucose dehydrogenase

Testing time of less than 10 seconds

Only a tiny blood sample of 0.6ul needed

200 test memory

7, 14, and 28 days averaging

Measurement range: 20-600 mg/dL

High-Low alert by colour

Strip ejector function in meter

Strip packed in individual foil

Long strip expiry of 24-months

Conforms to new ISO standard 15197:2013

Economical test strips

Large strip size

Custom designed lancing device

Blood sample: capillary whole blood

For in Vitro diagnostic use

For self-testing and professional use

Medical product

Size: 5.8 x 10.8 x 1.5 cm

Weight: 50g, including batteries

Batteries: 2 x AAA batteries

Battery life: 1000 tests

EAN no. : 747 2476 11000 3

Pack: 50/shipping carton. 5×10

