



PO-907

FINGER PULSE OXIMETER

Measurement of arterial oxygen saturation

(SpO2) and heart rate (Pulse)

For use at home and on the go

For persons with: Heart failure, Chronic bronchitis,
Asthma

Also for sportsman going at high altitude

OLED display

Low perfusion: 0.2%

Audio and video alarm function

Alarm and beep adjusted in settings menu

Pulse curve display

Automatic switch off

Anti-shaking

6 display modes

Big font display

Advanced DSP algorithm for high accuracy

Medical product

Incl. 2 x 1.5V AAA batteries

Carrying pouch

Size: approx. 3.4 x 6.2 x 3.1

Weight: approx. 50g, including batteries

Measuring range:

- SpO2: 70% - 100%

- Pulse: 30 - 250 beats/min

Batteries: 2 x AAA

EAN no.: 747 7869 90700 0

Pack: 100/Shipping carton

