

CERTEZA[®]
...your good health



FT-708

DIGITAL FLEXIBLE THERMOMETER

- Digital Flexible Thermometer
- Fast reading time
- Extra large display
- °F or °C settings
- Clinical accuracy
- Water proof
- Auto shut-off
- Fever Alert
- Break resistance
- Beeps when complete
- Memory recall
- For use at multiple sites
- BPA Free, Latex Free
- No mercury or glass
- Memory of last measured test
- With protective cover
- Medical product
- Size: approx. 2.2 x 13.9 x 1.3 cm
- Weight: approx. 11g, including battery
- Range: 90.0°F - 109.9°F
- Accuracy: $\pm 0.2^{\circ}\text{F}$ or $\pm 0.1^{\circ}\text{C}$
- Battery: 1 x 1.55V LR41
- EAN no.: 747 3830 70800 3
- Pack: 1 (10xFT-708)/ Shipping carton: 30x1

