

CERTEZA[®]
...your good health



BM-303

WRIST BLOOD PRESSURE MONITOR

Blood pressure and pulse measurement on wrist
Measurement during inflation
New attractive design
4 Colour backlit display
For wrist circumference of 13.5 - 19.5 cm
Movement and talk error indicator
Large LCD display
Hypertension risk indicator
2 x 199 memory recall function
Irregular Heartbeat Indicator
Average of last 3 readings
One touch automatic operation
Date and Time Function
Automatic switch off
Signal in case of application error
Low battery indicator
Storage case
Medical product

For wrist circumference: 13.5 - 19.5 cm
Size: approx. 8.4 x 6.6 x 2.6 cm
Weight: approx. 115g, excluding batteries
Measuring method: Oscillometric method
Measuring range:
- Blood pressure: 45 mmHg - 230mmHg
- Pulse: 40 - 199 beats/min
Batteries: 2 x AAA (included)
Case for storage

EAN no.: 747 9276 30300 0
Pack: 30 pcs / Shipping carton





BM-350

WRIST BLOOD PRESSURE MONITOR

Blood pressure and pulse measurement on wrist
Measurement during inflation
New attractive design
High definition display
Dual power source (AAA and Type C port)
Voice output of results
Dedicated toggle button for users
For wrist circumference of 13.5 - 19.5 cm
Movement and talk error indicator
Hypertension risk indicator
2 x 99 memory recall function
Irregular Heartbeat Indicator
Average of last 3 readings
One touch automatic operation
Date and Time Function
Automatic switch off
Signal in case of application error
Low battery indicator
Storage case
Medical product

For wrist circumference: 13.5 - 19.5 cm
Size: approx. 8.4 x 7.0 x 2.9 cm
Weight: approx. 127g, excluding batteries
Measuring method: Oscillometric method
Measuring range:
- Blood pressure: 45 mmHg - 230mmHg
- Pulse: 40 - 199 beats/min
Power Source: 2 x AAA (included) or Type C
Case for storage

EAN no.: 747 9276 35000 4
Pack: 30 pcs / Shipping carton



CERTEZA[®]
...your good health



BM-404

ARM BLOOD PRESSURE MONITOR

Blood pressure and pulse measurement on arm
Measurement during inflation
For arm circumference of 22-42 cm
USB-C external power source port
Movement error indicator
Cuff wrap OK alert
Large LCD display
Hypertension risk indicator
2 x99 memory recall function
Irregular Heartbeat Indicator
Average of last 3 readings
One touch automatic operation
Date and Time Function
Automatic switch off
Signal in case of application error
Low battery indicator
Storage pouch
Medical product

For arm circumference: 22 – 42 cm
Size: approx. 11.8 x 9.8 x 6.3 cm
Weight: approx. 220g, excluding batteries
Measuring method: Oscillometric method
Measuring range:
– Blood pressure: 57 mmHg – 255mmHg
– Pulse: 40 – 199 beats/min
Batteries: 4 x AAA (included)
USB-C adapter : 5V 1A (optional)
Pouch for storage

EAN no.: 747 8276 40400 7
Pack: 24 pcs / Shipping carton



CERTEZA[®]
...your good health



BM-470

ARM BLOOD PRESSURE MONITOR

Blood pressure and pulse measurement on arm
Measurement during inflation
Rechargeable Lithium battery
For arm circumference of 22-42 cm
USB-C charging port (cable included)
Movement error indicator
Cuff wrap OK alert
Voice output of results
Large LCD display
Hypertension risk indicator
2 x 120 memory recall function
Irregular Heartbeat Indicator
Average of last 3 readings
One touch automatic operation
Date and Time Function
Automatic switch off
Signal in case of application error
Low battery indicator
Storage Bag
Medical product

For arm circumference: 22 – 42 cm
Size: approx. 12.3 x 5.9 x 2.8 cm
Weight: approx. 225g
Measuring method: Oscillometric method
Measuring range:
– Blood pressure: 57 mmHg – 255mmHg
– Pulse: 40 – 199 beats/min
Power: 3.7 V Rechargeable Lithium battery
USB-C Type C cable included
Storage Bag

EAN no.: 747 8276 47000 2
Pack: 20 pcs / Shipping carton

