

NEW

Type C
Cuff Wrap alert
Movement indicator

CERTEZA[®]
...your good health



BM-407

ARM BLOOD PRESSURE MONITOR

Blood pressure and pulse measurement on arm
For arm circumference of 22-36 cm
Large LCD display
USB-C external power source port
Movement error indicator
Intelligent cuff wrap indicator
WHO classification of hypertension indicator
Irregular Heartbeat Indicator
Average of last 3 readings
One touch automatic operation
2×60 memory spaces
Medical product
Date and Time Function
Automatic switch off
Signal in case of application error
Low battery indicator
Zipper Storage case
AC adapter is optional

For arm circumference: 22-36 cm
Size: approx. 11.4 x 16.6 x 7.2 cm
Weight: approx. 425g, excluding batteries
Measuring method: Oscillometric method
Measuring range:
– Blood pressure: 0mmHg – 300mmHg
– Pulse: 30 – 180 beats/min
Batteries: 4 x AA

EAN no.: 747 8276 40700 8
Pack: 12/shipping carton

