

With
Length
Measurement
Scale

CERTEZA[®]
...your good health



BS-820

ELECTRONIC BABY SCALE

Baby scale with large side wall weighing platform
Length Measurement Scale
Large, easy to read white backlight LCD-display
Digit height 28.5 mm
Automatic weight lock function for moving babies
Tare weight function
5 g graduation (0.011 lb)
20 kg capacity (44 lb)
Automatic switch off
Easy to clean
Change over unit to kg/lb
Overload indicator
Batteries included (2 x AAA)
Gift box with practical handle

Size: approx. 563 x 333 x 88 mm
Range: Upto 20kg
Graduation: 5 grams
Accuracy: ± 0.011 lbs or ± 5 g
Digit size: 28.5mm
Battery: 2 x AAA

EAN no.: 747 3270 82000 6
Pack: 3 pcs / Shipping carton

