



**HP-220**

## HEATING PAD (PVC)

Heating pad with digital controller  
Consistent heat for optimal therapy  
Adjustable heat settings to control your therapy  
High-end controller with safety system  
9 Temperature settings  
3 Timer Settings (30min, 45min, 60min)  
Turbo heat up  
Button switch with LED lights  
30 x 40 cm detachable design for easy washing  
Soft, micro fleece material  
Automatic switch-off

Size: 30 x 40 cm  
Material: Micro fleece  
Power supply: 220V - 2240V 50Hz/ 60Hz  
Rated power: 52W  
Working temperature: 45-70°C

EAN no.: 747 4873 22000 2  
Pack: 12 pcs / Shipping carton

Limited warranty

