



BM-408

ARM BLOOD PRESSURE MONITOR

Blood pressure and pulse measurement on arm

Measurement during inflation

New attractive design

For arm circumference of 22-42 cm

USB-C external power source port

Movement error indicator

Intelligent cuff wrap indicator

Large LCD display

Hypertension risk indicator

2 x150 memory recall function

Irregular Heartbeat Indicator

Average of last 3 readings

One touch automatic operation

Date and Time Function

Automatic switch off

Signal in case of application error

Low battery indicator

Storage pouch

Medical product

For arm circumference: 22 – 42 cm

Size: approx. 14.98 x 8.04 x 4.57 cm

Weight: approx. 155g, excluding batteries

Measuring method: Oscillometric method

Measuring range:

– Blood pressure: 60 mmHg – 260mmHg

Diastolic 30-200mg

– Pulse: 30 – 180 beats/min

Batteries: 3 x AAA (included)

USB-C adapter : 5V 1A (optional)

Pouch for storage

EAN no.: 747 8276 40800 5

Pack: 12 pcs / Shipping carton

