



HP-250

HEATING PAD

Heating pad with digital controller
Consistent heat for optimal therapy
Adjustable heat settings to control your therapy
High-end controller with safety system
8 Temperature settings
3 Timer Settings (30min, 45min, 60min)
Turbo heat up
Button switch with LED lights
30 x 60 cm detachable design for easy washing
Soft, micro fleece material
Automatic switch-off

Size: 30 x 60 cm

Material: Micro fleece

Power supply: 220V - 2240V 50Hz/ 60Hz

Rated power: 63W

Working temperature: 35-65°C

EAN no.: 747 4873 25000 9

Pack: 20 pcs / Shipping carton

Suitable for all body areas like : Back, Thigh, Calf etc

*Limited warranty

