

CERTEZA®
...your good health



NB-608

NEBULIZER COMPRESSOR SYSTEM

Compressed air nebulization
Oil-free piston type
Storage space for nebulization kit
For treatment of upper and lower airways, colds,
asthma and respiratory diseases
Optimal respirable dose delivery
Quick effective treatment
Support for wide voltage input
Durable quality product
Non-slip body material
Built-in medicine cup holder
Angled-mouth piece
Power cord with fuse
Long tube
Medical product

Size: approx. 28.5 x 18.0 x 11.6 cm
Gross Weight: approx. 2.4 kg
Nebulizer chamber capacity: 6ml max
Nebulization rate: 0.27 ml/min
Particle size: 4 μ m (MMD) | 1.78 μ m (MMAD)
Fine Particle Fraction (%) (<4.7 μ m) : 92 %
Noise level: Around 57 dBA
Electrical requirement: 230V / 50Hz
Power consumption: Below 65 watts
Max. Compressor Pressure: 35 to 50 psi
Liter Flow Range: 5-8 LPM

EAN no.: 747 6322 60800 3
Pack: 6/shipping carton

